

# Medlink CDR

Compact, Advanced X-Ray System Designed for  
Chiropractic Care That Includes New 45° Arm Rotation



**MEDLINK** IMAGING

a **VIEWWORKS** company


# Compact Chiropractic X-Ray System with 45° Arm Rotation

## Medlink CDR, Advanced X-Ray Technology

Medlink has a new upgraded compact x-ray system designed specifically for chiropractic care. The Medlink CDR has an ultra small footprint and comes with a variety of features such as full software integration, **45° arm rotation**, magnetic lock brakes, 40kW generator, dual LED laser collimator and a removable grid.

This x-ray solution can fit in virtually any x-ray room as well as drastically cut down installation time reducing clinical downtime.

**1-Year Warranty**  
(extended warranty available)



**The VXVUE acquisition software** is an award-winning intelligent software solution which utilizes the most advanced imaging processing. VXvue's image viewer program contains specialized modules to acquire, adjust, and manage diagnostic images created by the VIVIX-S Series.



**Bucky Enclosure**  
Features fully counterbalanced powerful magnetic lock brakes from the handle control to raise and lower the system effortlessly. The bucky also rotates to maintain tube alignment during arm angulation positioning.

**MEDLINK** I M A G I N G

---

a viewworks company